

bCONTM Biocontainment System



Two bCONTM Racks . . . Four bCONTM Cages . . . Maximum Species Flexibility

In the research community, requirements of the modern Biocontainment facility constantly evolve in response to new knowledge and opportunity. For decades, the most widely accepted animal caging used in ABSL-3 facilities has been the simple, static (non-ventilated) Micro-IsolatorTM—and for good reason . . . until we introduced the **bCONTM** Biocontainment System. Combining strict engineering controls, flexibility and superior performance with ease of use, **bCONTM** provides a simple and reliable ABSL-3 animal housing solution that dramatically improves upon static caging performance.

And now, you can provide a complete Multi-speciesTM high-density housing for all your rodents with only two racks, the **bCONTM** Super Mouse 1800TM and **bCONTM** One Cage 2100TM housing systems. For **mice**, select the Super Mouse 1800TM and mix and match Super Mouse 750TM Micro-IsolatorTM cages in any combination you choose. For maximum flexibility, you can house **mice, rats, hamsters** or **guinea pigs** with interchangeable One Cage 2100TM and One CageTM Micro-IsolatorTM cages on one rack, featuring modular snap-on feeders.

- Negative cages in controlled environment cabinet for housing mice in ABSL-3
- Dual-Layer Protection: mechanical & air curtain barriers between cages & room
- Containment & Micro-Isolation: cage to cage, cage to room & room to cage
- AirlockTM Cage Exhaust with negative pressure—cages do not become positive
- Cage filter top Locking Security with quick & quiet operation
- Cage Docking Security with automatic visual indicator
- Micro-ShieldTM Cage Filter prolongs HEPA Filter
- Available UPS battery power system maintains rack ventilation during power outages
- During airflow disruption cages revert to static mode—Animals can breathe
- Minimal components for ease of use & excellent cage visibility
- Decontaminate entire system: Chlorine Dioxide Gas or Vaporized Hydrogen Peroxide
- Maximum Species Flexibility with a choice of two racks and four cages
- **bCONTM** Super Mouse 1800TM Interchangeable Micro-IsolatorTM Housing System
Super Mouse 750TM **bCONTM** Micro-IsolatorTM (Floor Area >75 in² / >483 cm²)
Super Mouse 1800TM **bCONTM** Micro-IsolatorTM (Floor Area >180 in² / >1161 cm²)
- **bCONTM** One Cage 2100TM Interchangeable Micro-IsolatorTM Housing System
One Cage 2100TM **bCONTM** Micro-IsolatorTM (Floor Area >210 in² / >1354 cm²)
One CageTM **bCONTM** Micro-IsolatorTM (Floor Area >80 in² / >516 cm²)



*b*CON™ Biocontainment System

Effective Dual-Layer Protection for Your People and Your Animals

Specifically engineered for maximum species flexibility in ABSL-3 facilities, the *b*CON™ Biocontainment System design places a negative cage within a controlled environment cabinet system providing dual-layer protection—at the cage and cabinet level. Engineering controls create a 3-zone containment system, the *b*CON™ negative Micro-Isolator™ Cage (primary), the controlled environment of the cabinet (secondary) and the animal holding room (tertiary).


Cage-Level Containment & Cage-to-Cage Airflow Balance

The *b*CON™ Biocontainment System maintains Micro-Isolator™ cages exclusively under negative pressure; surrounded by HEPA filtered cabinet air, cages do not become positive. HEPA filtered cabinet air is drawn through each cage filter top membrane (indirect connection) replacing the same amount of air evacuated through our unique Airlock™ cage exhaust connection (direct connection). With any number of docked cages on the rack, the system provides cage-to-cage airflow balance—without compromising performance—with simple, effective cage-level containment, previously unachieved with ABSL-3 animal housing.

While docking or removing cages, the *b*CON™ Biocontainment System airflow patterns continue to protect the integrity of the Micro-Isolator™ environment. Insert a cage into the fully docked position and a hinged plate in front of the cage moves downward for visual confirmation. During power outages, an optional UPS back-up battery power system maintains rack ventilation. In a worse case scenario, cages revert to a static Micro-Isolator™ condition; for decades, the most widely accepted type of animal housing used in ABSL-3.




One Cage™



One Cage 2100™



Super Mouse 750™



Super Mouse 1800™

E²™ Environmental Control Systems

E²™ Cabinet Supply and E²™ Cage Exhaust Units provide microprocessor controlled HEPA filtered (99.99% efficient to .3 micron particles) airflow, protecting the *b*CON™ Biocontainment System and are interfaced through a touch screen control panel with digital menus. Operational parameters are continuously monitored including; cage air changes, cabinet airflow, air pressure, temperature, humidity, fan performance and HEPA filter hours. Remote monitoring, historical reporting and alarm alerts can be accessed with a facility network or up to 60 days of data stored into on-board memory is available for download onto a memory card.

