

# New RAIR Enviro-Gard™ Cage Monitor

Featuring self-calibrating ventilated rack technology.

## New & Improved.

- Self-Calibrating, Plug & Play Rack Certification
- Validates Cage-To-Cage Airflow Balance (ACH)
- Compensates for Altitude Location Pressure Variance
- Monitors and Displays ACH, Temperature, Time & Date
- Eliminates Need for Rack Certification Service Call
- Data Logging & USB Download of Spreadsheet Data
- Rechargeable Battery or Power Adapter Operation
- Easy Integration with the Wireless Alert System
- Available for Super Mouse 750™ Super Mouse 1800™, World Cage 500™, One Cage™, One Cage 2100™, and Super Rat 1400™ Micro-Isolator® Systems.



Monitor, record and transfer key micro-environmental data, with the reliability you need to ensure your animal housing units are running smoothly. Certify cage-level airflow balance with self-calibrating rack technology in just minutes per rack, with consistent measurements you can trust. You can even self-calibrate your racks while they are populated with rodents.

The New Cage Monitor provides a precision monitor/certification instrument that is a cost-effective, easy-to-configure, low maintenance solution that makes your work simpler. Record and monitor Air Changes per Hour (ACH), Temperature, Time and Date from any ventilated Micro-Isolator® cage location. Data can be stored onboard the cage unit, downloaded directly to your PC (with Alert System integration) or USB drive and monitored remotely in real-time with an internet or network connection.

Now you can know precisely what's going on inside your AllerZone™ Micro-Isolator® cages with reliable technology built upon decades of biomedical research equipment design experience.

trusted engineering. better research. **lab products**