

Bench-Mate™

Mini Workstation

Designed to assist with your day-to-day laboratory work, Bench-Mate™ provides an invaluable aid for your animal manipulations and cage changing procedures. A retractable platform extends over mouse cages, bringing your work closer while improving posture. Three separate disinfection wells locate your instruments in the center of operations, reducing time and motion. This design allows convenient disinfection of the platform, forceps and other instruments between aseptic procedures.



- Reduces time and motion
- Retractable platform
- Improves posture
- Separate disinfection wells
- Compact, portable design
- Reduces bench-top clutter
- Vacuum grip base

Bench-Mate™ is the perfect ergonomic companion for enhancing your animal manipulations and cage changing procedures. The slender 5.25" wide by 8.50" high profile conserves bench-top space and reduces bench-top clutter. The retractable 6.75" by 6.50" stainless steel grid platform provides a handy staging surface. And the powerful vacuum grip base enables you to relocate Bench-Mate™ at your convenience and secure it to most any flat surface, even wet surfaces.

Benchmate™

lab products

