

Enviro-Gard™ A Environmental Control System

With the Enviro-Gard™ A System, balanced environmental control and monitoring of HEPA filtered airflow is easier than ever.

Enviro-Gard™ A combines microprocessor and volumetric airflow sensor technology to control and maintain precise airflow. Even as the HEPA filter becomes loaded, fan speed is automatically adjusted to compensate. Consistent airflow is provided for each ventilated Micro-Isolator™ rack, regardless of the number of cages being housed. Each unit is equipped with a large LED digital display of airflow, LCD digital display with scrolling menu of various operational parameters, green and red indicator lights of key operational status, and much more. At your option, Enviro-Gard™ A can be networked to provide remote monitoring capabilities and automatic email or pager alerts.



Enviro-Gard™ A Air Supply and Exhaust Units

Compact and lightweight, energy efficient with low noise and low vibration, Enviro-Gard™ A delivers state of the art environmental controls designed to take the guesswork out of equipment operation and performance.

- 24 Hour airflow control
- Continuously monitors airflow, filters, fan, room temperature
- LED digital display of airflow (ACH/CFM)
- LCD digital display of operational parameters
- LED status indicators of filters and fan
- Energy efficient with low noise and vibration
- Supports email and pager warnings
- Interfaces with Enviro-Gard™ i Software
- UL, CSA and CE certification



The
Enviro-Gard™
System



Enviro-Gard™ B Environmental Control System

Enviro-Gard™ B combines microprocessor and volumetric airflow sensor technology to control and maintain precise airflow. Even as the HEPA filter becomes loaded, fan speed is automatically adjusted to compensate. Consistent airflow is provided for each ventilated Micro-Isolator™ rack, regardless of the number of cages being housed. Each unit is equipped with an LED digital display of airflow. HEPA filter, pre-filter and fan status is indicated by green and red indicator lights.

Enviro-Gard™ B Air Supply and Exhaust Units

Compact and light weight, energy efficient with low noise and low vibration, Enviro-Gard™ B delivers state of the art environmental controls designed to take the guesswork out of equipment operation and performance.

- 24 Hour airflow
- Continuously monitors airflow, filters fan, room temperature
- LED digital display of airflow (ACH/CFM)
- LED status indicators of filters and fan
- Energy efficient with low noise and vibration
- UL, CSA and CE certification



The
Enviro-Gard™
System

