

# Enviro-Gard™ MultiPlex™

## Environmental Control Systems

### Multiple Rack, HEPA Filtered Air Handling Unit



One MultiPlex™ handles up to 10 ventilated racks at once!

- Maintains Precise Airflow with Microprocessor Controls
- Takes the Guesswork Out of Operation and Performance
- Handles Different Styles of Ventilated Racks / Cages Simultaneously
- Capacity: 400 CFM Output (May Vary with Selection of Duct Design)
- Monitors Airflow, Filters and Fan
- Large LED Digital Display of Accurate Airflow (ACH & CFM or CMH)
- LCD Digital Display of numerous Operational Parameters
- LED Status (Red / Green Indicators) of HEPA filter, Pre-filter and Fan
- Environmental Monitoring Option - Enviro-Gard™ i Facility Integration Software
- Supports Remote Monitoring, Reporting, Email Text Alerts
- Configurations: Supply Only, Exhaust Only or Supply & Exhaust
- Available in Free Standing, Wall Mount or Mobile Designs
- Exhaust Connection Option with Facility HVAC Systems
- Large HEPA Filter for Extended Life
- Energy Efficient with Low Noise and Low Vibration

The  
*Enviro-Gard*<sup>TM</sup>  
System