


Modular Primate Caging Unit

INFINITE POSSIBILITIES!

Standard Features & Options

Choose from this partial listing of features and options or contact us to discover more.



Conveniently attach units together to transfer your animals horizontally for transport or preparing for cage wash processing. The one-over-one design provides full vertical access for animal exercise or group housing flexibility. Each cage compartment features a single pull, suspended squeeze back with a double bar perch and optional Delrin platform pad. Self storing, sliding socialization panels provide animals with; Separation, Viewing, Finger Touching or Full Access with a solid panel rear refuge. Meet or exceeds current USDA requirements for housing Group 4 Primates (10-15 kg).

- **Interlocking Unit Socialization:**
 - Interlock (1/1) Unit to Unit continuously
 - Multiple Cage Configurations limited only by number of interlocked units per row
 - 4 Stage Socialization, Adjacent Cages
 - Expansion Horizontal and / or Vertical
- **Doors:**
 - Full Swing, Removable Cage Front with Guillotine Door (std.)
 - Optional (within Full Swing Front):
 - Sliding Door,
 - Dual Sliding Doors
 - Specialty Doors
- **Pans:**
 - Dry Pan (SST or Plastic)
 - Optional:
 - Rear Flush (to wall, trough)
 - Rear Flush w/Downspout
 - Metabolism Collection
- **Body Panels:**
 - Clear Polycarbonate Top & Rear (std.)
 - Interchangeable Mesh Side (std.)
 - Optional:
 - SST Rear Panels
 - Mesh Top Panel
- **Feeders & Watering:**
 - Box Feeder (std.)
 - Optional:
 - Foraging (Reach Through Turf),
 - Puzzle / Foraging (Top Mounted)
 - Bottle Watering,
 - Automatic Watering
- **Floors:**
 - Two Removable Grid Floors (std.)
 - Remove Mid Section Floor, Solid Panel and Pan for Full Vertical Unit Expansion
- **Squeeze Back:**
 - Single Pull/Perch design featuring two perch tubes and a SST platform
 - Optional: Perch Platform Pad
- **Enrichment: (Adjacent Cage)**
 - 4-Stage Socialization Panels (std.);
 - Separation: Self-Storing, Solid, Sliding
 - Viewing: Self-Storing, Sliding w/View Panel
 - Finger Touching: Dual Mesh Fixed Panels
 - Full Access: Horiz. Cage Access & Retreat

