

One Cage 2100™ Stack-A-Rack™

Maximum Flexibility



Stack-A-Rack™ Interchangeable See-Through™ Housing Systems meet or exceed ILAR Guide Recommendations for housing rodents.

Create your own stackable housing system with our unique Stack-A-Rack™ Interchangeable See-Through™ Housing Systems. Housing rats, mice, hamsters and guinea pigs, the One Cage 2100™ Stack-A-Rack™ provides the ultimate in flexibility and convenience by housing two distinctly different cage sizes. Simply rotate cage runner supports to switch from the larger and smaller cages. House both simultaneously on the same rack. Up to seven modular tiers accept as many as 42 One Cage 2100™ or 84 One Cage™ cages—or anywhere in between!

Stack-A-Rack™ systems provide a variety of ordering and configuration options. To meet various housing density and installation requirements, single-sided and double-sided mobile systems are available. You can order Stack-A-Rack™ unassembled, preassembled or with welded, fixed shelving. One Cage 2100™ and One Cage™ cages accept LifeSpan™ rodent enrichment, providing a stimulating multilevel cage environment. Both cages also accept Micro-Isolator™ Filter Tops, allowing cage inventory to be shared with animals housed in Micro-Isolator™ Systems, for other housing system applications.



One Cage 2100™
(210 in² floor area)



One Cage™
(80 in² floor area)

- Mount 1 to 7 modular shelves to create your own flexible system
- Rack includes a rigid mobile base, angle uprights and 1 to 7 modular shelves
- Components can be disassembled and cleaned / autoclaved separately
- Two Species-flexible™ interchangeable cages
- One Cage 2100™ floor area >210 in² (>1354 cm²)
- One Cage™ floor area >80 in² (>516 cm²)
- Cages share interchangeable, snap-on feeders to adapt for various species
- Eliminates wire bar lids for easy observations of animals, food and water
- Faster cage changes with minimal components and reduced inventory
- Watering options: Hydropac™, water bottles or automatic watering

