

One Cage™

AllerZone™ Micro-Isolator™

High Density Ventilated Racks

Independently
Tested
& Certified



Universal Housing For All Your Small Rodents.

- Maximizes Capacity (up to 112 Cages)
- One Cage™ Floor Area ($>80 \text{ in}^2 / >516 \text{ cm}^2$)
- One Flexible System for Mice, Rats, Hamsters and Guinea Pigs
- Adapts for Various Species with Modular, Snap-On Feeders
- AllerZone™ Control of Airborne Allergens
- Positive Cage / Negative Rack Protection of Animals and Personnel
- Precise Airflow Balance from Cage to Cage (+/- .1 CFM)
- Airflow Remains Balanced when Rack is Partially Full of Cages
- Micro-Isolator™ Effectiveness is Not Compromised if Power Loss Occurs
- Adequate Airflow for Animals Remains if Power Loss Occurs
- Eliminates Wire Bar Lids for Easy Observations of Animals, Food and Water
- Faster Cage Changes with Direct Access to Animals and Reduced Inventory
- Watering Options: (**New**) Hydropac™, Bottles or Automatic
- (**New**) LifeSpan™ Rodent Enrichment
- Cage Monitor Control Unit Option - Cage Level Monitoring / Control / Certification
- Environmental Monitoring Option - Enviro-Gard™ i Facility Integration Software
- Connects with Facility HVAC Exhaust and / or Supply Systems
- Meets or Exceeds ILAR Guide Recommendations

Lab Products, Inc. • 742 Sussex Avenue, P.O. Box 639, Seaford, DE 19973
800.526.0469 • 302.628.4300 • Fax 302.628.4309 • labproductsinc.com

© Bio Medic Corporation 2007. One Cage, AllerZone, Micro-Isolator, LifeSpan and Enviro-Gard are trademarks of Lab Products, Inc.
Hydropac is a trademark of Hydropac / Lab Products, Inc. Covered by various U.S. and International Patents and Patents pending. Made in USA.