

Stay-Clean™ Duplex Workbench

What's new in Class 100 Workbenches?

The Stay-Clean™ Duplex Workbench provides a Class 100, HEPA filtered (99.99% efficient to 0.3 micron particles) environment for the manipulation, handling and transferring of animals to clean cages. HEPA filtered air flows down from the top of the workbench, pulling airborne particulate from the cages and into the air intake around the perimeter of the work surface. This airflow reduces personnel exposure to particulate escaping from the cages and helps protect animals from personnel.

This mobile unit features an easy to view digital display, providing pre-filter and fan status. The unit is compact, may be operated from either side, and has an adjustable work surface height for ergonomic flexibility.

- Compact and mobile
- Double sided
- Streamlines workflow
- Digital controls
- Adjustable work surface height 32"-41"

*The Stay-Clean™ Duplex Workbench . . .
providing the right solutions and the right
environment for researchers and the
animals they care for!*



Appropriate aseptic techniques are required to help assure protection of the animals. Not approved for use with biohazard, toxic or flammable materials.

