

Stay-Clean™ Workbenches

Laminar Flow Work Environments



First and Foremost, Designed To Meet Your Changing Needs.

- Maintains Precise Airflow with Digital Microprocessor Controls
- Takes the Guesswork Out of Operation and Performance
- HEPA Filtered Air / 99.99% efficient to .3 Micron Particles
- Large Color Coded LED Digital Display of HEPA Filter Status
- LED Status (Red / Green Indicators) of Pre-filter and Fan Status
- Ergonomic, Push Button Height Adjustment (ELD Models)
- Adjustable Work Surface Height from 37" to 46" (93.98 to 116.84 cm)
- Adjustable Sash Opening from 12" to 17" (30.48 to 43.18 cm)
- Helps Protect Animals from Cross Contamination
- Helps Minimize Personnel Exposure to Particulate
- Connects with Facility HVAC Exhaust Systems (8" Collar Adapter Req'd.)
- Overall Dimensions (Digital Controls):
#30909-D 59.13" x 34" x 76.50" H (150.2 x 86.4 x 194.3 cm H)
#30910-D 83.13" x 34" x 76.50" H (211.2 x 86.4 x 194.3 cm H)
- Overall Dimensions (Digital Controls and Adjustable Height):
#30909-ELD 59.13" x 34" x 76.5"/ 85.5" H (150.2 x 86.4 x 194.3 / 217.2 cm H)
#30910-ELD 83.13" x 34" x 76.5"/ 85.5" H (211.2 x 86.4 x 194.3 / 217.2 cm H)