

Stay-Clean™ Workbenches.

Laminair Flow Work Environments.
First and Foremost, Designed To Meet Your Changing Needs.

- Maintains Precise Airflow with Digital Microprocessor Controls
- Takes the Guesswork Out of Operation and Performance
- HEPA Filtered Air / 99.99% efficient to .3 micron particles
- Large Color Coded LED Digital Display of HEPA Filter Status
- LED Status (red/green indicators) of Pre-filter and Fan Status
- Ergonomic, Push Button Height Adjustment (ELD Models)
- Adjustable Work Surface Height from 37" to 46" (93.98 to 116.84 cm)
- Adjustable Sash Opening from 12" to 17" (30.48 to 43.18 cm)
- Helps Protect Animals From Cross Contamination
- Helps Minimize Personnel Exposure to Particulate
- Connects with Facility HVAC Exhaust Systems (8" Collar Adapter Required.)



30910-ELD

Stay-Clean™ Workbench Mobile Work Station With Digital Display & Electric Lift

Workbench Catalog Number	Class 100 Airflow	Overall Dimensions (in. / mm)			Work Surface Dimensions (in. / mm)				Pre-filter	HEPA Filter	Exhaust HEPA Filter	² Air Velocity	³ Power Requirements
		Width	Depth	¹ Height	Surface Area / Adjustable Height								
					Width	Depth	Min.	Max.					
30909-ELD	Horizontal Positive	59.13 1501.90	34.00 863.60	76.50 1943.10	49.38 1254.25	25.25 641.35	37.00 939.80	46.00 1168.40	30977-B 2 per unit	30982	30982-EX	100 (±10) LFPM	110 VAC 50/60 Hz 15 Amps
30910-ELD	Horizontal Positive	83.13 2111.50	34.00 863.60	76.50 1943.10	73.38 1863.85	25.25 641.35	37.00 939.80	46.00 1168.40	30978-B 3 per unit	30983	30983-EX		

¹85.50" / 2171.70 mm at maximum lift height. ²LFPM (Linear Feet per Minute). ³Can be factory set for 220 VAC 50/60 Hz operation. Sash opening adjusts from 12.00" / 304.80 mm (down) to 16.50" / 419.10 mm (up).



30909-ELD



30909-D



30910-D

Stay-Clean™ Workbench Mobile Work Station With Digital Display

Workbench Catalog Number	Class 100 Airflow	Overall Dimensions (in. / mm)			Work Surface (in. / mm)			Pre-filter	HEPA Filter	Exhaust HEPA Filter	¹ Air Velocity	² Power Requirements
		Width	Depth	Height	Surface Area		Surface Height					
					Width	Depth						
30909-D	Horizontal Positive	59.13 1501.90	34.00 863.60	76.50 1943.10	49.38 1254.25	25.25 641.35	37.00 939.80	30977-B 2 per unit	30982	30982-EX	100 (±10) LFPM	110 VAC 50/60 Hz 15 Amps
30910-D	Horizontal Positive	83.13 2111.50	34.00 863.60	76.50 1943.10	73.38 1863.85	25.25 641.35	37.00 939.80	30978-B 3 per unit	30983	30983-EX		

¹LFPM (Linear Feet per Minute). ²Can be factory set for 220 VAC 50/60 Hz operation. Sash opening adjusts from 12.00" / 304.80 mm (down) to 16.50" / 419.10 mm (up).