

Super Mouse 750™

Designed exclusively for housing mice.

**Independently
Tested
& Certified**

The Super Mouse 750™ High Density Housing System now holds up to 900 adult mice with 180 cages per rack, in virtually the same footprint! Conceived to enhance your productivity, each component provides time saving benefits beginning with the "no wire bar lid" design, allowing faster cage changes and easier observations of the animals, food and water at a glance. The clear plastic Modular Diet Delivery System accepts either a Hydropac™ water pouch or a water bottle. As an option, your research studies can benefit from LifeSpan™ rodent enrichment while keeping your animals happier, healthier and drier.



AllerZone™ technology is ***Independently Tested and Certified**. Super Mouse 750™ includes AllerZone™ technology, creating a separate pressure zone between the cage and the room, while preventing cage air from escaping to the room and room air from entering the cage. The positive pressure cage protects animals with HEPA filtered air, and airflow of the negative pressure rack directs animal dander, particulate, and odors away from personnel, providing allergen control. AllerZone™ is independently tested for containment of Mouse Urine Protein (M.U.P.) and SF₆ Tracer Gas.

Standard systems are available to house from 63 up to 180 cages per rack and are fully compatible with previously designed Micro-Isolator LP™ and Super Mouse Micro-Isolator™ cages.

- Maximizes housing capacity
- Accepts Hydropac™ water pouches
- Eliminates old style wire bar lids
- Allows faster cage changes
- Provides easier observations
- Precise airflow balance
- Helps control airborne allergens
- Meets ILAR Guide recommendations



lab products



*The Efficiency of Ventilated Cage Racks as Engineering Controls in Reducing Personnel Exposure to Mouse Urine Protein. John Schafer, MFS, CIH, HEIM, CEPA. Gas Challenge Testing performed by NSF accredited technicians of ENV Services, Inc., Testing and Certification of contamination Control Equipment.

Lab Products, Inc. • 742 Sussex Avenue, P.O. Box 639, Seaford, DE 19973
800.526.0469 • 302.628.4300 • Fax 302.628.4309 • labproductsinc.com

© Bio Medic Corporation 2008. Super Mouse 750, LifeSpan, AllerZone, Micro-Isolator LP and Micro-Isolator are trademarks of Lab Products, Inc. Hydropac is a trademark of Hydropac / Lab Products, Inc. Covered by various U.S. and International Patents and Patents pending. Made in USA.

