


Super Rat 1400™

AllerZone™
Micro-Isolator™ System



Universal Housing For All Your Small Rodents.

- Super Rat 1400™ Floor Area (>140 in² / >903 cm²)
- Effective in Ventilated (Rack) and Non-Ventilated (Ambient) Conditions
- AllerZone™ Control of Airborne Allergens (Ventilated)
- One Flexible System for Rats, Mice, Hamsters and Guinea Pigs
- Eliminates Wire Bar Lids for Easy Observations of Animals, Food and Water
- Faster Cage Changes with Direct Access to Animals and Reduced Inventory
- Snap-Together Filter Top Design for Easy Filter Changing
- Modular Diet Delivery System with Snap-On Feeder Attachments
- Modular Diet Delivery System can be Positioned at Cage Rear or Cage Front
- Designed to Maximize Internal Cubic Area Available for Animals
- Cage Divider Option Provides Separated / Dual Compartments
- Watering Options: (**New**) Hydropac™, Bottles or Automatic
- Available in Several Materials Including Polycarbonate:
 - Zyfone™ (2 year warranty)
 - Zytem™ (5 year warranty)
- Meets or Exceeds ILAR Guide Recommendations

Lab Products, Inc. • 742 Sussex Avenue, P.O. Box 639, Seaford, DE 19973
800.526.0469 • 302.628.4300 • Fax 302.628.4309 • labproductsinc.com

© Bio Medic Corporation 2007. Super Rat 1400, AllerZone, Micro-Isolator, Zyfone and Zytem are trademarks of Lab Products, Inc.
Hydropac is a trademark of Hydropac / Lab Products, Inc. Covered by various U.S. and International Patents and Patents pending. Made in USA.