


Wireless Alert System

Stay informed about your ventilated racks, cages and research animals around the clock, wirelessly.

Peace of Mind.

- Alert Monitoring Around the Clock
- Instant Service & Alarm Messages
- Easy to Install, Use and Maintain
- Remote Access to Data - Anytime
- In-Depth Real-Time System Analysis
- Download & Store Spreadsheet Data
- Wireless Savings Over Cable System
- Reduces Physical Inspection Time



Works with any device with an internet connection, including Smartphones, PDA's, Tablets, Laptops and more!



The Alert System allows you to remotely monitor the supply and exhaust HEPA filtered air systems of your ventilated racks, while sending electronic alerts to personnel in the event of a fault condition. Messages can be sent by e-mail or cell phone text, even during a power outage.

Login to view data anytime through your web-enabled computer or hand held device. Download data for spreadsheet archiving and regulatory agency reports. Once logged in, the status of all connected racks will be displayed with an intuitive, interactive screen. An optional wall mounted display panel can be mounted in a hallway for monitoring a suite of rooms, with highly visible red/green LED operation status indicators and interactive screen.

Stay in touch with what's going on inside your Micro-Isolator® cages by linking with the optional Cage Monitor for recording and monitoring Air Changes per Hour (ACH), Temperature, Time and Date from any ventilated Micro-Isolator® cage location.

trusted engineering. better research. **lab products**