

# Super Mouse 750™ Micro-Isolator® Systems



## *Designed Exclusively for Housing Your Laboratory Mice*

Provide high density housing for your mouse colonies with the Super Mouse 750™ AllerZone™ Micro-Isolator® System... designed with minimal components to simplify handling, reduce inventory and increase your productivity.

- Floor Area: >75 in<sup>2</sup> / >483 cm<sup>2</sup>
- Effective in Ventilated (Rack) & Non-Ventilated (Ambient) Conditions
- Patented AllerZone™ Control of Airborne Allergens
- No "Wire Bar Lid" - Clear Observations, Easier Cage Changes & Handling
- Easy Filter Changing with a Snap-Together Filter Top
- LifeSpan™ Rodent Enrichment adds >18 in<sup>2</sup> surface roaming area
- Interchangeable with \*Econo-Cage® on Your Ventilated Racks
- Modular Diet System can be Positioned at Cage Rear or Front
- Cage Divider Option Provides Separated, Dual Compartments
- Designed to Maximize Internal Cubic Area Available for Rodents
- Watering: **(NEW)** Hydropac®, Bottles or Automatic
- Available in Several Materials Including Polycarbonate; Zyfone™ (2 year warranty) & Zytem™ (5 year warranty)
- Meets or Exceeds ILAR Guide Recommendations
- Many More Options Available - Contact Lab Products for Details



*LifeSpan™ Rodent Enrichment and other components are interchangeable with all Super Mouse 750™ and Super Mouse 1800™ Micro-Isolator® Systems.*

\*Econo-Cage® Disposable Micro-Isolator® System - contact us for details.

**lab products**

Lab Products, Inc. • 742 Sussex Avenue, P.O. Box 639, Seaford, DE 19973 • 800.526.0469 • 302.628.4300 • Fax 302.628.4309 • labproductsinc.com

© Bio Medic Corporation 2012. Super Mouse 750, Super Mouse 1800, AllerZone, LifeSpan, Zyfone and Zytem are trademarks of Lab Products, Inc. Micro-Isolator, Econo-Cage and Hydropac are registered trademarks of Lab Products, Inc., Bio Medic Corporation and Hydropac / Lab Products, Inc. respectively. Covered by various U.S. and International Patents and Patents pending. Made in USA.



# Super Mouse 750™ Micro-Isolator® Systems

(CAGE FLOOR AREA: >75 in <sup>2</sup> / >483 cm <sup>2</sup> )							ORDERED SEPARATELY				
MATERIAL	MICRO-ISOLATOR® SYSTEM CATALOG NO.			CAGE BOTTOM	ALLERZONE™ FILTER TOP	MODULAR DIET DELIVERY SYSTEM	LIFESPAN™ RODENT ENRICHMENT	WATER BOTTLE 15 oz. 443 ml.	ONE CAP™	SIPPER TUBE	NEW HYDROPAC® WATERING
	STATIC	VENTILATED									
	NO GROMMET	AIR GROMMET	AIR & WATER GROMMET								
Polycarb.	75050PC	-GA	-GAW	75001	75051PC	75006-SL	75009PC	75040	30211	30110S	Contact Lab Products for details.
Zyfone™	75050ZF	-GA	-GAW	75031	75051ZF	75034-SL	75009ZF	75043	30211ZF	30110S	
Zytem™	75050ZT	-GA	-GAW	75021	75051ZT	75024-SL	<sup>1</sup> 75009ZF	75042	30211ZT	30110S	

1) LifeSpan™ is available in Polycarbonate or Zyfone™ materials.

Mice			
>25 g	up to 25 g	up to 15 g	<10 g
5	6	9	12



75051PC

75009PC



LifeSpan™ Rodent Enrichment



75054 Filter Sheet



75006-SL



75051PC-GA



75050PC-GA



75057 Super Mouse 750™  
Filter Replacement Device



30509H  
Zytem™ Cardholder  
Horizontal Flip-Up



30508  
SS Cardholder  
Horizontal  
Clip-On, Flip-Up



30508V  
SS Cardholder  
Vertical  
Clip-On, Flip-Up



Hydropac®  
Water Pouch



30110S  
30211  
75040  
Water Bottle



-GAW  
Air & Water Grommet

lab products

Lab Products, Inc. • 742 Sussex Avenue, P.O. Box 639, Seaford, DE 19973 • 800.526.0469 • 302.628.4300 • Fax 302.628.4309 • labproductsinc.com

© Bio Medic Corporation 2012. Super Mouse 750, Super Mouse 1800, AllerZone, LifeSpan, Zyfone and Zytem are trademarks of Lab Products, Inc. Micro-Isolator, Econo-Cage and Hydropac are registered trademarks of Lab Products, Inc., Bio Medic Corporation and Hydropac / Lab Products, Inc. respectively. Covered by various U.S. and International Patents and Patents pending. Made in USA.