

rabbit metal cage units



40364



40317

## Removable cages suspended over individual excreta trays

Product Number	Description	Overall Dimensions W x D x H Inches	Cage Size		Animals per Cage	Food Hopper	Watering	Excreta Tray	Casters <sup>(1)</sup> Inches
			Dimensions W x D x H Inches	Usable Floor Area Sq. Feet					
40317 <sup>(2)</sup>	6 removable cages suspended on 3-tier rack (Wire-sided cages) <sup>(3)</sup>	57 x 28 x 70½	26 x 26 x 17¼	4	One up to 5.4 kg.	"J" Type	Bottle, can or automatic <sup>(4)</sup>	Tray rests on rack under each cage	5
40318	8 removable cages suspended on 4-tier rack (Wire-sided cages) <sup>(3)</sup>	57 x 28 x 80	26 x 26 x 15¼	"	"	"	"	"	"
40319	8 removable cages suspended on 4-tier rack (Wire-sided cages) <sup>(3)</sup>	49 x 28 x 80	22 x 26 x 15¼	3	One up to 4 kg.	"	"	"	"
40320	9 removable cages suspended on 3-tier rack (Wire-sided cages) <sup>(3)</sup>	72½ x 28 x 70½	22 x 26 x 17¼	"	"	"	"	"	"
40363	6 removable cages suspended on 3-tier rack	56 x 32 x 70½	24¼ x 30¼ x 17	5	One over 5.4 kg.	"	"	"	"
40364	6 removable cages suspended on 3-tier rack	56 x 26 x 70½	24¼ x 24¼ x 17	4	One up to 5.4 kg.	"	"	"	"
27130	4 removable cages suspended on 2-tier rack	58 x 28 x 68¼	26 x 26 x 16¼	"	"	"	"	Metabolism tray	"
25951	6 removable cages suspended on 3-tier rack	56½ x 27¼ x 70½	24½ x 25¼ x 17¼	"	"	"	"	Tray rests on rack under each cage	"
31610	8 removable cages suspended on 4-tier rack	59 x 36 x 85	27 x 33 x 15½	5	One over 5.4 kg.	"	"	"	"

(1) For caster information, see pages 63 and 115.

(2) Available with modifications for housing poultry.

(3) Available also with perforated metal sides. For details, contact Lab Products.

(4) For automatic watering systems, contact Lab Products.

Note: To order replacement components, contact Lab Products.

# Cages and racks integrated

Product Number	Description	Overall Dimensions W x D x H Inches	Cage Size				Animals Accommodated		Food Hopper	Watering	Excreta Tray	Casters <sup>(1)</sup> Inches
			Individual (divider in)		Tier (divider out)		Individual Cage	Tier				
			Dimensions W x D x H Inches	Usable Floor Area Sq. Feet	Dimensions W x D x H Inches	Usable Floor Area Sq. Feet						
40305	4-Tier rack with 8 integrated cages	36 x 30 x 78¼	17 x 29¾ x 15¼	3	34 x 29¾ x 15¼	6	1, up to 4 kg.	4, up to 2 kg. 2, 2-4 kg. 1, over 4 kg.	"J" Type	Bottle or automatic <sup>(2)</sup>	One tray rests on rack under 2 cages	5
40324	"	48 x 30 x 78¼	23 x 29¾ x 15¼	4	46 x 29¾ x 15¼	8	1, up to 5.4 kg.	5, up to 2 kg. 2, 2-5.4 kg. 1, over 5.4 kg.	"	"	"	"
40330	"	42 x 30 x 78¼	20 x 29¾ x 15¼	3.5	40 x 29¾ x 15¼	7	1, up to 4 kg.	4, up to 2 kg. 2, 2-4 kg. 1, over 4 kg.	"	"	"	"
40331	"	38 x 30 x 77¼	17½ x 27½ x 14	3	35 x 27½ x 14	6	"	"	"	"	"	8

(1) For caster information, see pages 63 and 115.

(2) For automatic watering systems, contact Lab Products.

Note: To order replacement components, contact Lab Products.



## Removable cages totally self-contained

Product Number	Description	Overall Dimensions W x D x H Inches	Cage Size		Animals per Cage	Food Hopper	Watering	Excreta Tray	Casters <sup>(1)</sup> Inches
			Dimensions W x D x H Inches	Usable Floor Area Sq. Feet					
40316 <sup>(2)</sup>	6 totally self-contained cages, suspended on 3-tier rack	57 x 24 x 76½	24¼ x 24 x 20½	4	One up to 5.4 kg.	"J" Type	Bottle or automatic <sup>(3)</sup>	Integrated into cage	5
40307	As above, cages resting on rack	63 x 28 x 73½	30 x 28 x 19¾	5	One over 5.4 kg.	"	"	"	"

(1) For caster information see pages 63 and 115.  
 (2) Available with modifications for housing poultry.  
 (3) For automatic watering systems, contact Lab Products.

## Individual cage

Product Number	Description	Cage Size		Animals per Cage	Food Hopper	Watering	Excreta Tray	Rack
		Dimensions W x D x H Inches	Usable Floor Area Sq. Feet					
41302	Totally self-contained cage	24 x 28 x 18	4	One up to 5.4 kg.	"J" Type	Bottle	Integrated into cage	See page 118



40307



40316