



Econo-Cage®

DISPOSABLE SYSTEM



Mouse Cage

- Low Cost
- Easy to Use
- Stores Compactly
- Eliminates Washing
- BPA Free (Bisphenol A)
- Fits Super Mouse Racks
- Ventilated or Static
- Recyclable



Econo-Cage® System
with Air Grommet.

*The task of housing rodents is complicated enough.
The Econo-Cage® system simplifies the entire process.*



Uses **Hydropac®**, Water Bottles
or Automatic Watering!

Only Two Parts

From a truly exceptional design, two things stand out ... the cage top and the cage bottom. There's nothing more to handle.

Simple Is Better

A seamlessly integrated diet unit and cage top provides excellent visibility and reduced overall handling.

Smooth Edges - Safer Handling

Made from precision injection molded FDA Grade polypropylene plastic, the Econo-Cage® is available standard or irradiated.

Featuring AllerZone™ Technology

The Econo-Cage® is completely compatible with Super Mouse Ventilated Racks & Cages, and it's the only disposable cage with AllerZone™ technology. If you already use Super Mouse housing, you don't have to buy new racks to use the Econo-Cage®.

lab products inc

Econo-Cage® DISPOSABLE SYSTEM



72123

Econo-Cage® Disposable System

- 72123 Econo-Cage® Disposable System
- 72123-B Econo-Cage® Disposable System w/ 1/4" Deep Bedding
- 72123-I Econo-Cage® Disposable System (Irradiated)
- 72328-BI Econo-Cage® Disposable System w/ 1/4" Deep Bedding (Irradiated)



72456-G

Econo-Cage® Disposable System w/ Air Grommet

- 72456-G Econo-Cage® Disposable System
- 72456-GB Econo-Cage® Disposable System w/ 1/4" Deep Bedding
- 72456-GI Econo-Cage® Disposable System (Irradiated)
- 72456-GBI Econo-Cage® Disposable System w/ 1/4" Deep Bedding (Irradiated)



72456-FG

Econo-Cage® Disposable System w/ Flap Grommet

- 72456-FG Econo-Cage® Disposable System
- 72456-FGB Econo-Cage® Disposable System w/ 1/4" Deep Bedding
- 72456-FGI Econo-Cage® Disposable System (Irradiated)
- 72456-FGBI Econo-Cage® Disposable System w/ 1/4" Deep Bedding (Irradiated)



72001

Econo-Cage® Disposable Cage Bottom

Shown w/ 1/4" Bedding

- 72001 Econo-Cage® Disposable Cage Bottom
- 72001-B Econo-Cage® Disposable Cage Bottom w/ 1/4" Deep Bedding
- 72001-I Econo-Cage® Disposable Cage Bottom (Irradiated)
- 72001-BI Econo-Cage® Disposable Cage Bottom w/ 1/4" Deep Bedding (Irradiated)



72014-G

Econo-Cage® Disposable Cage Bottom w/ Air Grommet

Shown w/ 1/4" Bedding

- 72014-G Econo-Cage® Disposable Cage Bottom
- 72014-GB Econo-Cage® Disposable Cage Bottom w/ 1/4" Deep Bedding
- 72014-GI Econo-Cage® Disposable Cage Bottom (Irradiated)
- 72014-GBI Econo-Cage® Disposable Cage Bottom w/ 1/4" Deep Bedding (Irradiated)



72014-FG

Econo-Cage® Disposable Cage Bottom w/ Flap Grommet

Shown w/ 1/4" Bedding

- 72014-FG Econo-Cage® Disposable Cage Bottom
- 72014-FGB Econo-Cage® Disposable Cage Bottom w/ 1/4" Deep Bedding
- 72014-FGI Econo-Cage® Disposable Cage Bottom (Irradiated)
- 72014-FGBI Econo-Cage® Disposable Cage Bottom w/ 1/4" Deep Bedding (Irradiated)

Available in the following quantities: Case of 100 / Pallet of 1,600 / Lot of 12,800.

72020

Econo-Cage® Reusable Cardholder



- 72020 Econo-Cage® Reusable Cardholder

Available in a Case of 100.

lab products